

TOOL CRIB STATION



FEATURES

- Automate your existing tool room, crib, or store
- Issue/return tools without an attendant
- Capture transaction before unlocking door to crib
- Use at remote locations in a "runner system"

PHYSICAL

- Touch screen or keypad
- Handheld high-output barcode laser scanner
- Full one-year warranty

SOFTWARE

- Locate by item or location: easily search for any item in any tool crib in the AutoCrib system, using AutoCrib Software
- Perform physical or cycle counts
- Purchasing and receiving: perform stocking transactions
- All transactions recorded in real-time
- Windows OS
- EDI capable (thru AutoCrib software)
- Tool Crib Station uses AutoCrib Software (web-based AutoCrib.net or client/server software). The software features an Intelligent Dashboard with key metrics visible at a glance on one screen. Over 200 reports are pre-loaded.

You can quickly and easily convert your manual tool crib to a virtual dispensing machine, with AutoCrib Tool Crib Station. Just designate an existing tool room, tool crib or stores area as a tool crib, using AutoCrib Software. Install an AutoCrib Tool Crib Station at the crib window — and now, tool crib attendants or users perform the data entry using a touch screen. After the user selects the desired items, a receipt with bin locations is printed to assist the attendant or user in locating the items within the secured area.

You can also place Tool Crib Stations at various locations or remote work cells throughout a plant to create a "runner system". Users enter tool or supplies requests at the Tool Crib Station. The central tool crib prints out the order, and issues and delivers the requested tool by a runner (at lower cost), while keeping the operator at his work center.

HOW IT WORKS

Select the desired function on the touch screen or keypad (Issue, Return, Locate, etc.). To issue an item from the tool crib, enter an employee number and pin number. The user can also be prompted to enter cost center data, such as department or job number, if applicable.

Select the item desired. The User can enter the part number or name of the item desired, and then select the item from a list. A photo can also be viewed to confirm the correct item has been selected. The user selects the quantity desired. At this point, the crib door can be unlocked by the system. A receipt with bin numbers is printed to assist in locating the items.

To return an item, enter an employee number and pin number. A list of items to be returned by that employee appears. Select the item to be returned. All transactions are recorded in the AutoCrib database in real-time

Runner System also available